
commentrgb0.5,0.0,0.0 keywordrgb0.0,0.5,0.0 keywordtypergb0.38,0.25,0.125
keywordflowrgb0.88,0.5,0.0 preprocessorrrgb0.5,0.38,0.125 stringliteralrgb0.0,0.125,0.25
charliteralrgb0.0,0.5,0.5 vhldigitrgb1.0,0.0,1.0 vhdlkeywordrgb0.43,0.0,0.43 vhdl-
logicrgb1.0,0.0,0.0 vhdlcharrgb0.0,0.0,0.0
darkgray

analyzer

0.1

Generated by Doxygen 1.8.7

Sat Mar 11 2017 12:18:15

Contents

Chapter 1

Module Index

1.1 Modules

Here is a list of all modules:

Main menu	??
Calibration mode	??
Etalon's data write	??

Chapter 2

Class Index

2.1 Class List

Here are the classes, structs, unions and interfaces with brief descriptions:

current_t	??
etalon_t	??

Chapter 3

File Index

3.1 File List

Here is a list of all files with brief descriptions:

/home/zvebabi/Documents/workspace/LEDMicrosensor_project0/analizer/ analizer.cpp	??
/home/zvebabi/Documents/workspace/LEDMicrosensor_project0/analizer/ analizer.h	??
/home/zvebabi/Documents/workspace/LEDMicrosensor_project0/analizer/Release/ analizer.d	??
/home/zvebabi/Documents/workspace/LEDMicrosensor_project0/analizer/Release/ SPI.d	??

Chapter 4

Module Documentation

4.1 Main menu

Functions

- void `loop` ()

4.1.1 Detailed Description

4.1.2 Function Documentation

4.1.2.1 void `loop` ()

Infinite cycle Will perform after startup device Available in main menu commands:

- **q**
Go to factory calibration function.
- **e**
Express calibration method. Not implemented yet
- **s**
One mode of measurements with option of coeffs calibration
Manually choosing preamps parameters
Sample&Holde mode
- **m**
One mode of measurements with option of coeffs calibration
Manually choosing preamps parameters
Inpulse measurements mode
- **a**
One mode of measurements with option of coeffs calibration
Series of measurements with predefined preamp parameters
Sample&Holde mode

Definition at line 81 of file analyzer.cpp.

4.2 Calibration mode

Functions

- void `factoryCalibr()`

4.2.1 Detailed Description

4.2.2 Function Documentation

4.2.2.1 void `factoryCalibr()`

Main calibration function Available in calibration menu commands:

- **w**
Set width of impulse in microseconds
step is 1us, max 200us
- **t**
Save to eeprom value of input resistance of second OpAmp
- **a**
Set currents for LED in channel 1 Max value is 1000mA, with step 1mA
- **b**
Set currents for LED in channel 2 Max value is 1000mA, with step 1mA
- **g**
reset LEDs if something wrong
- **f**
Light-up first led
- **n**
Next led choose
- **o**
Start infinite series of pulse to LED
- **p**
Make one pulse to LED
- **i**
Disable LED
- **z**
Return current led number
- **c**
Set coeffs of OpAmp
Control digital potentiometers ad5141
Max value is 100kOhm, with step 0.39kOhm

- **r**
Go to preAmp calibration
- **s**
Save parameters for all LED to EEPROM
- **e**
Exit from calibration mode

Definition at line 149 of file analyzer.cpp.

4.3 Etalon's data write

Functions

- void `preAmpCalibr` ()

4.3.1 Detailed Description

4.3.2 Function Documentation

4.3.2.1 void `preAmpCalibr` ()

Amplifier calibration function Available in PreAmp calibration menu commands:

- **c**
Set coeffs of OpAmp
Control digital potentiometers ad5141
Max value is 100kOhm, with step 0.39kOhm
- **x**
Choose etalon cell in memory
- **e**
Exit from calibration mode
- **s**
Save parameters to EEPROM
- **z**
Return current etalon cell number
- **k**
Save k for etalon
- **m**
set min and max value of g2
- **n**
Compute Ai
- **a**
Compute Ai with doMeasurementsSH_Avg method
- **l**
Write Ai in manual mode

Definition at line 374 of file analyzer.cpp.

Chapter 5

Class Documentation

5.1 `current_t` Struct Reference

```
#include <analyzer.h>
```

Public Attributes

- `uint16_t curr1`
- `uint16_t curr2`

5.1.1 Detailed Description

Definition at line 65 of file analyzer.h.

5.1.2 Member Data Documentation

5.1.2.1 `uint16_t current_t::curr1`

Definition at line 67 of file analyzer.h.

5.1.2.2 `uint16_t current_t::curr2`

Definition at line 68 of file analyzer.h.

The documentation for this struct was generated from the following file:

- `/home/zvebabi/Documents/workspace/LEDMicrosensor_project0/analyzer/analyzer.h`

5.2 `etalon_t` Struct Reference

```
#include <analyzer.h>
```

Public Attributes

- `float k [NUM_OF_LED]`
- `float g1`

- float [g2mid](#)
- float [g2min](#)
- float [g2max](#)

5.2.1 Detailed Description

Definition at line 53 of file analyzer.h.

5.2.2 Member Data Documentation

5.2.2.1 float etalon_t::g1

Definition at line 56 of file analyzer.h.

5.2.2.2 float etalon_t::g2max

Definition at line 59 of file analyzer.h.

5.2.2.3 float etalon_t::g2mid

Definition at line 57 of file analyzer.h.

5.2.2.4 float etalon_t::g2min

Definition at line 58 of file analyzer.h.

5.2.2.5 float etalon_t::k[**NUM_OF_LED**]

Definition at line 55 of file analyzer.h.

The documentation for this struct was generated from the following file:

- [/home/zvebab/ Documents/ workspace/ LEDMicrosensor_ project0/ analyzer/ analyzer.h](#)

Chapter 6

File Documentation

6.1 /home/zvebabi/Documents/workspace/LEDMicrosensor_project0/analizer/analizer.cpp File Reference

```
#include "analizer.h"
```

Functions

- [ISR](#) (TIMER1_COMPA_vect)
- int [main](#) (void)
- void [setup](#) ()
- void [loop](#) ()
- void [factoryCalibr](#) ()
- void [preAmpCalibr](#) ()
- void [doMeasurements](#) (uint8_t numOfEtalon, bool calcNorm)
- void [doMeasurementsSH](#) (uint8_t numOfEtalon, bool calcNorm)
- void [doMeasurementsSH_Avg](#) (bool calcNorm)
- void [readADCOneTime](#) (uint16_t &value)
- void [readADC](#) (float &value)
- void [writeConfigToUart](#) ()
- void [setC_R](#) (float val)
- void [doOnePulse](#) (uint16_t pulseWidth)
- void [dischargeSampleHold](#) ()
- void [shiftRegisterReset](#) ()
- void [shiftRegisterNext](#) ()
- void [shiftRegisterFirst](#) ()
- void [disableLED](#) ()
- void [SerialClean](#) ()
- void [setPreAmp](#) (float RWB1, float RWB2)
- void [setCurrent](#) (uint8_t channelN, uint16_t curValue)
- void [setPulseWidth](#) (uint16_t width)

6.1.1 Function Documentation

6.1.1.1 void disableLED ()

Utility function for calibration mode

Definition at line 983 of file analizer.cpp.

6.1.1.2 void dischargeSampleHold ()

Discharge capasitor, while measurements going in Sample&Hold mode

Definition at line 923 of file analyzer.cpp.

6.1.1.3 void doMeasurements (uint8_t numOfEtalon = 0, bool calcNorm = false)

One mode of measurements with option of coeffs calibration Manually choosing preamps parameters Inpulse measurements

Parameters

|>p0.10|>p0.15|p0.678|

in *numOfEtalon* - which etalon we use for fast calibration

in *calcNorm* - true - calibration, false - measurements

Definition at line 568 of file analyzer.cpp.

6.1.1.4 void doMeasurementsSH (uint8_t numOfEtalon = 0, bool calcNorm = false)

One mode of measurements with option of coeffs calibration Manually choosing preamps parameters Sample&Holde mode

Parameters

|>p0.15|p0.805|

numOfEtalon - which etalon we use for fast calibration

calcNorm - true - calibration, false - measurements

Definition at line 648 of file analyzer.cpp.

6.1.1.5 void doMeasurementsSH_Avg (bool calcNorm = false)

One mode of measurements with option of coeffs calibration Series of measurements with predefined preamp parameters Sample&Holde mode

Parameters

|>p0.15|p0.805|

calcNorm - true - calibration, false - measurements

Definition at line 728 of file analyzer.cpp.

6.1.1.6 void doOnePulse (uint16_t pulseWidth)

Make one pulse to the LED with previously setted current and time of pulse

Definition at line 915 of file analyzer.cpp.

6.1.1.7 ISR (TIMER1_COMPA_vect)

Analyzer.cpp

Author

LED Microsensor

Definition at line 8 of file analizer.cpp.

6.1.1.8 int main (void)

Definition at line 15 of file analizer.cpp.

6.1.1.9 void readADC (float & value)

Utility function for Sample&Holde mode of measurements make 16 samples and calc average,

Parameters

|>p0.10|>p0.15|p0.678|

out *value* buffer for reading voltage (mV)

Definition at line 823 of file analizer.cpp.

6.1.1.10 void readADCOneTime (uint16_t & value) [inline]

Utility function for inpulse mode of measurements

Parameters

|>p0.10|>p0.15|p0.678|

out *value* buffer for reading adc value (w/o convert to mV)

Definition at line 814 of file analizer.cpp.

6.1.1.11 void SerialClean ()

Utility function for cleanserial port buffer

Definition at line 991 of file analizer.cpp.

6.1.1.12 void setC_R (float val)

Save to eeprom value of input resistanse of second OpAmp

Parameters

|>p0.15|p0.805|

val

Definition at line 910 of file analizer.cpp.

6.1.1.13 void setCurrent (uint8_t channelN, uint16_t currValue)

This function set currents for LED in channel channelN Control DAC Max value is 1000mA, with step 1mA After current will set, it send a message with currently setted values to serial port

Parameters

|>p0.15|p0.805|

channelIN - Number of channel 1 or 2

currValue - DAC output voltage, max 1000mV it is equal 1000mA on LED

Definition at line 1044 of file analyzer.cpp.

6.1.1.14 void setPreAmp (float *RWB1*, float *RWB2*)

This function set coeffs of OpAmp Control digital potentiometers ad5141 Max value is 100kOhm, with step 0.39kOhm After resistance will set, it send a message with currently setted values to serial port

Parameters

|>p0.15|p0.805|

RWB1 1st cascade (current to voltage)

RWB2 2nd cascade (signal amplification)

Definition at line 1004 of file analyzer.cpp.

6.1.1.15 void setPulseWidth (uint16_t *width*)

Set width of impulse in microseconds step is 1us, max 200us

Parameters

|>p0.15|p0.805|

width - time of pulse, max 200us

Definition at line 1089 of file analyzer.cpp.

6.1.1.16 void setup ()

Main initialization of MCU

Definition at line 27 of file analyzer.cpp.

6.1.1.17 void shiftRegisterFirst ()

Utility function for work with shift register Select first pair of led

Definition at line 970 of file analyzer.cpp.

6.1.1.18 void shiftRegisterNext ()

Utility function for work with shift register Select next led to ON

Definition at line 959 of file analyzer.cpp.

6.1.1.19 void shiftRegisterReset ()

Utility function for work with shift register Set all outputs to High level

Definition at line 934 of file analyzer.cpp.

6.1.1.20 void writeConfigToUart ()

Print all saved to eeprom data to serial port

Definition at line 843 of file analizer.cpp.

6.2 /home/zvebabi/Documents/workspace/LEDMicrosensor_project0/analizer/analizer.h File Reference

```

#include "Arduino.h"
#include <stdio.h>
#include <util/delay.h>
#include <avr/eeprom.h>
#include <SPI.h>

```

Classes

- struct [etalon_t](#)
- struct [current_t](#)

Macros

- #define [GeneratorPin](#) 9
- #define [ShiftRegisterDelay](#) 1
- #define [DACDelay](#) 1
- #define [WBDelay](#) 1000
- #define [DISCHARGE_DELAY](#) 50
- #define [PULSE_DELAY](#) 25
- #define [REFERENCE_V](#) 3000.0
- #define [NUM_OF_LED](#) 45
- #define [NUM_OF_ETALON](#) 3
- #define [MAX_PULSE_WIDTH](#) 150
- #define [MAX_CURRENT](#) 1000
- #define [ADC_PORT](#) PORTE
- #define [DAC_PORT](#) PORTF
- #define [PREAMP_PORT](#) PORTB
- #define [SS_ADC](#) PE6
- #define [SS_PREAMP](#) PB4
- #define [SS_DAC](#) PF1
- #define [SR_ENABLE](#) PB7
- #define [SR_CLR](#) PD4
- #define [SR_DATA](#) PD7
- #define [SR_CLK](#) PD6
- #define [SH_PORT](#) PORTD
- #define [SH_SET](#) PD0
- #define [SH_RESET](#) PD1
- #define [GEN1](#) OCR1A
- #define [GEN2](#) OCR1B

Functions

- void [setPreAmp](#) (float RWB1, float RWB2)
- void [setCurrent](#) (uint8_t channelIN, uint16_t currValue)
- void [setPulseWidth](#) (uint16_t width)
- void [setC_R](#) (float val)
- void [doOnePulse](#) (uint16_t pulseWidth)
- void [dischargeSampleHold](#) ()
- void [disableLED](#) ()
- void [factoryCalibr](#) ()
- void [preAmpCalibr](#) ()
- void [doMeasurementsSH](#) (uint8_t numOfEtalon=0, bool calcNorm=false)
- void [doMeasurementsSH_Avg](#) (bool calcNorm=false)
- void [doMeasurements](#) (uint8_t numOfEtalon=0, bool calcNorm=false)
- void [readADCOneTime](#) (uint16_t &value)
- void [readADC](#) (float &value)
- void [writeConfigToUart](#) ()
- void [shiftRegisterReset](#) ()
- void [shiftRegisterNext](#) ()
- void [shiftRegisterFirst](#) ()
- void [SerialClean](#) ()
- void [setup](#) ()
- void [loop](#) ()
- void [__cxa_pure_virtual](#) ()

Variables

- uint8_t EEMEM [_empty](#) [20] = {0xF}
- uint16_t EEMEM [_pulseWidth](#) = 80
- float EEMEM [_c_R](#) = 3.9
- [current_t](#) EEMEM [_pairsOfCurrent](#) [NUM_OF_LED]
- [etalon_t](#) EEMEM [_etalons](#) [NUM_OF_ETALON]
- float EEMEM [_coefficients](#) [NUM_OF_LED]
- volatile [current_t](#) [cur4AllLed](#) [NUM_OF_LED]
- volatile float [coeffs](#) [NUM_OF_LED]
- volatile bool [oneTimes](#) = false
- volatile bool [pulseEnd](#) =false

6.2.1 Macro Definition Documentation

6.2.1.1 `#define ADC_PORT PORTE`

Definition at line 28 of file analyzer.h.

6.2.1.2 `#define DAC_PORT PORTF`

Definition at line 29 of file analyzer.h.

6.2.1.3 `#define DACDelay 1`

Definition at line 17 of file analyzer.h.

6.2.1.4 **#define DISCHARGE_DELAY 50**

Definition at line 19 of file analizer.h.

6.2.1.5 **#define GEN1 OCR1A**

Definition at line 47 of file analizer.h.

6.2.1.6 **#define GEN2 OCR1B**

Definition at line 48 of file analizer.h.

6.2.1.7 **#define GeneratorPin 9**

Analizer.cpp

Author

LED Microsensor

Definition at line 15 of file analizer.h.

6.2.1.8 **#define MAX_CURRENT 1000**

Definition at line 25 of file analizer.h.

6.2.1.9 **#define MAX_PULSE_WIDTH 150**

Definition at line 24 of file analizer.h.

6.2.1.10 **#define NUM_OF_ETALON 3**

Definition at line 23 of file analizer.h.

6.2.1.11 **#define NUM_OF_LED 45**

Definition at line 22 of file analizer.h.

6.2.1.12 **#define PREAMP_PORT PORTB**

Definition at line 30 of file analizer.h.

6.2.1.13 **#define PULSE_DELAY 25**

Definition at line 20 of file analizer.h.

6.2.1.14 **#define REFERENCE_V 3000.0**

Definition at line 21 of file analizer.h.

6.2.1.15 #define SH_PORT PORTD

Definition at line 42 of file analyzer.h.

6.2.1.16 #define SH_RESET PD1

Definition at line 44 of file analyzer.h.

6.2.1.17 #define SH_SET PD0

Definition at line 43 of file analyzer.h.

6.2.1.18 #define ShiftRegisterDelay 1

Definition at line 16 of file analyzer.h.

6.2.1.19 #define SR_CLK PD6

Definition at line 39 of file analyzer.h.

6.2.1.20 #define SR_CLR PD4

Definition at line 37 of file analyzer.h.

6.2.1.21 #define SR_DATA PD7

Definition at line 38 of file analyzer.h.

6.2.1.22 #define SR_ENABLE PB7

Definition at line 36 of file analyzer.h.

6.2.1.23 #define SS_ADC PE6

Definition at line 31 of file analyzer.h.

6.2.1.24 #define SS_DAC PF1

Definition at line 33 of file analyzer.h.

6.2.1.25 #define SS_PREAMP PB4

Definition at line 32 of file analyzer.h.

6.2.1.26 #define WBDelay 1000

Definition at line 18 of file analyzer.h.

6.2.2 Function Documentation

6.2.2.1 void __cxa_pure_virtual ()

Definition at line 224 of file analizer.h.

6.2.2.2 void disableLED ()

Utility function for calibration mode

Definition at line 983 of file analizer.cpp.

6.2.2.3 void dischargeSampleHold ()

Discharge capasitor, while measurements going in Sample&Hold mode

Definition at line 923 of file analizer.cpp.

6.2.2.4 void doMeasurements (uint8_t numOfEtalon = 0, bool calcNorm = false)

One mode of measurements with option of coeffs calibration Manually choosing preamps parameters Inpulse measurements

Parameters

|>p0.10|>p0.15|p0.678|

in *numOfEtalon* - which etalon we use for fast calibration

in *calcNorm* - true - calibration, false - measurements

Definition at line 568 of file analizer.cpp.

6.2.2.5 void doMeasurementsSH (uint8_t numOfEtalon = 0, bool calcNorm = false)

One mode of measurements with option of coeffs calibration Manually choosing preamps parameters Sample&Holde mode

Parameters

|>p0.15|p0.805|

numOfEtalon - which etalon we use for fast calibration

calcNorm - true - calibration, false - measurements

Definition at line 648 of file analizer.cpp.

6.2.2.6 void doMeasurementsSH_Avg (bool calcNorm = false)

One mode of measurements with option of coeffs calibration Series of measurements with predefined preamp parameters Sample&Holde mode

Parameters

|>p0.15|p0.805|

calcNorm - true - calibration, false - measurements

Definition at line 728 of file analizer.cpp.

6.2.2.7 void doOnePulse (uint16_t *pulseWidth*)

Make one pulse to the LED with previously setted current and time of pulse

Definition at line 915 of file analyzer.cpp.

6.2.2.8 void readADC (float & *value*)

Utility function for Sample&Holde mode of measurements make 16 samples and calc average,

Parameters

|>p0.10|>p0.15|p0.678|

out *value* buffer for reading voltage (mV)

Definition at line 823 of file analyzer.cpp.

6.2.2.9 void readADCOneTime (uint16_t & *value*) [inline]

Utility function for inpulse mode of measurements

Parameters

|>p0.10|>p0.15|p0.678|

out *value* buffer for reading adc value (w/o convert to mV)

Definition at line 814 of file analyzer.cpp.

6.2.2.10 void SerialClean ()

Utility function for cleanserial port buffer

Definition at line 991 of file analyzer.cpp.

6.2.2.11 void setC_R (float *val*)

Save to eeprom value of input resistanse of second OpAmp

Parameters

|>p0.15|p0.805|

val

Definition at line 910 of file analyzer.cpp.

6.2.2.12 void setCurrent (uint8_t *channelN*, uint16_t *currValue*)

This function set currents for LED in channel *channelN* Control DAC Max value is 1000mA, with step 1mA After current will set, it send a message with currently setted values to serial port

Parameters

|>p0.15|p0.805|

channelN - Number of channel 1 or 2

currValue - DAC output voltage, max 1000mV it is equal 1000mA on LED

Definition at line 1044 of file analyzer.cpp.

6.2.2.13 void setPreAmp (float *RWB1*, float *RWB2*)

This function set coeffs of OpAmp Control digital potentiometers ad5141 Max value is 100kOhm, with step 0.39kOhm After resistance will set, it send a message with currently setted values to serial port

Parameters

|>p0.15|p0.805|

RWB1 1st cascade (current to voltage)

RWB2 2nd cascade (signal amplification)

Definition at line 1004 of file analizer.cpp.

6.2.2.14 void setPulseWidth (uint16_t *width*)

Set width of impulse in microseconds step is 1us, max 200us

Parameters

|>p0.15|p0.805|

width - time of pulse, max 200us

Definition at line 1089 of file analizer.cpp.

6.2.2.15 void setup ()

Main initialization of MCU

Definition at line 27 of file analizer.cpp.

6.2.2.16 void shiftRegisterFirst ()

Utility function for work with shift register Select first pair of led

Definition at line 970 of file analizer.cpp.

6.2.2.17 void shiftRegisterNext ()

Utility function for work with shift register Select next led to ON

Definition at line 959 of file analizer.cpp.

6.2.2.18 void shiftRegisterReset ()

Utility function for work with shift register Set all outputs to High level

Definition at line 934 of file analizer.cpp.

6.2.2.19 void writeConfigToUart ()

Print all saved to eeprom data to serial port

Definition at line 843 of file analizer.cpp.

6.2.3 Variable Documentation

6.2.3.1 float EEMEM _c_R = 3.9

Definition at line 76 of file analyzer.h.

6.2.3.2 float EEMEM _coefficients[NUM_OF_LED]

Definition at line 79 of file analyzer.h.

6.2.3.3 uint8_t EEMEM _empty[20] = {0xF}

allocate eeprom variable

Definition at line 74 of file analyzer.h.

6.2.3.4 etalon_t EEMEM _etalons[NUM_OF_ETALON]

Definition at line 78 of file analyzer.h.

6.2.3.5 current_t EEMEM _pairsOfCurrent[NUM_OF_LED]

Definition at line 77 of file analyzer.h.

6.2.3.6 uint16_t EEMEM _pulseWidth = 80

Definition at line 75 of file analyzer.h.

6.2.3.7 volatile float coeffs[NUM_OF_LED]

Definition at line 84 of file analyzer.h.

6.2.3.8 volatile current_t cur4AllLed[NUM_OF_LED]

Definition at line 83 of file analyzer.h.

6.2.3.9 volatile bool oneTimes = false

Definition at line 85 of file analyzer.h.

6.2.3.10 volatile bool pulseEnd =false

Definition at line 86 of file analyzer.h.

6.3 /home/zvebabi/Documents/workspace/LEDMicrosensor_project0/analizer/↵
Release/analizer.d File Reference

6.4 /home/zvebabi/Documents/workspace/LEDMicrosensor_project0/analizer/↵
Release/SPI.d File Reference